

Authorized Users

While an important part of your SOP, this list of all Authorized Users, including their badge number, should appear as Appendix C of this document and users should sign the document to verify that they have read and understand the Standard Operating Procedure for the LCA. Appendix C can be updated as frequently as needed.

Incidental Personnel

In addition to the authorized primary users, incidental personnel may be in the room at the time of the experiment. These personnel are not trained on the system, nor have they gone through the laser safety program. They will be observing only and not involved in the operation of the system. They will not be in the room during any alignment or maintenance of the lasers.

Normal Laser Operation

Room STRF 252, Greehey Campus contains the following lasers:

Class	Type	Make/Model	SN	Wavelength	Output Power mW
3B	Solid state	Coherent/1170506	L05060102	407 nm	100
3B	Solid state	Coherent/ Sapphire 488-200 CW	LDP.1137970.135015	488 nm	200
3B	Solid state	Coherent/ Sapphire 561-150 CW	LDP.1162896.132873	561 nm	150
3B	Fiber Laser	MPB/2RU-VFL-P-300-647	VFL-1104-1320	647 nm	300

All lasers are connected to the Nikon Ti microscope via a fiber optic cable. All lasers and their power supplies are combined into a wheeled sled that is placed on the floor behind the anti-vibration table. The sole purpose of the lasers is for super-resolution imaging.

Eyewear section

Wavelength specific eyewear will be used by field service engineers during alignment.

Alignment Hazard Control

All lasers are aligned by service engineers and are not adjusted by users or facility staff.

Laser Hazard Control

1. Access to laser control room is restricted to trained personnel. The room is secured by a card-access door reader. All access is approved by the campus police.
2. The entire system is enclosed by laser proof curtain.
3. Operation of all Class 3B lasers is only accessible through the system software. Computer access is restricted to trained users with unique login names and passwords.
4. A "Laser in Use" warning sign is mounted above the entrance of the enclosed system booth.
5. All lasers are enclosed with no open beam throughout the path.
6. All users are trained in the operation of the NSTORM and the proper use and care for the lasers.

Primary users are trained in the operation of the NSTORM and the proper use care for the lasers.

Control of Additional LCA hazards

Indicate other possible hazards associated with the lasers in your LCA.

Associated Chemical Hazard Control

List chemicals used in this LCA include a list MSDS numbers, or attach MSDS's to the end of the document (not required to attach MSDS sheets). If you prefer, provide the chemical list section from your Project Review Document as an Appendix. Indicate in this section if there is a registered Satellite Waste Accumulation Area in the LCA and where it is located. Discuss any site-specific chemical hazards for this LCA in this section.

Emergency Procedures

1. Shut down the laser system. Use "emergency stop" button if equipped.
2. Provide for the safety of the personnel, i.e. first aid, CPR, etc.
3. If a fire has been created as a result of the laser, follow appropriate procedures to put out the fire:
 - Within the surgical field, douse with sterile water
 - For fires in other areas, utilize an appropriate class fire extinguisher to extinguish the fire
4. Obtain medical assistance. In the event of a suspected eye injury, have the injured person keep their head upright and still to restrict any bleeding in the eye.
 - For life-threatening injuries (major burns, cardiac arrest following electrocution), dial 911 for immediate medical assistance.
 - For non-life threatening injuries (laser eye injuries, minor skin burns), employees should be evaluated by a physician as soon as possible. **Do not allow anyone with a potential laser eye injury to drive themselves.**

This SOP has been adapted with permission from Argonne National Laboratories.

User Name

Aaron Gonzalez
Breeanne Soteros
Brian Stoveken
Hui Yan
Kimi Kong
Lauren Schneider
Suman Kanji

User Login Email

gonzaleza33@livemail.uthscsa.edu
soteros@livemail.uthscsa.edu
stoveken@uthscsa.edu
yanh3@uthscsa.edu
kongk@uthscsa.edu
schneider1@uthscsa.edu
kanji@uthscsa.edu