

University of Texas Health Science Center at San Antonio  
Standard Operating Procedures For

Sweptfield Confocal  
Optical Imaging Core Facility  
STRF 252, Greehey Campus  
Laser Controlled Area Laser Controlled Area  
ver 4.0 (October 14, 2021)

Prepared by: Exing Wang

Laser Custodian:

Exing Wang \_\_\_\_\_  
Name Printed Signature Date

Laser Safety Officer:

**Jennifer Cerecero** \_\_\_\_\_  
Name Printed Signature Date

## INTRODUCTION

The Sweptfield confocal contains the following lasers:

Class	Type
3B	Solid state (405 nm)
3B	Solid state (488 nm)
3B	Solid state (561nm)
3B	Solid state (640 nm)

All lasers and their power supplies are integrated into a wheeled unit that is placed on the floor under the anti-vibration table. The output is connected to the Nikon Ti microscope via a fiber optic cable.

### Laser Custodian

Exing Wang                      53532                      (210)562-4062                      wange3@uthscsa.edu

-----  
Name                      Badge#                      Contact#                      Email

### Alternative Contact:

Jimmy Wewer                      38506                      (210)567-3151                      wewer@uthscsa.edu  
James Lechleiter                      28297                      (210)562-4043                      lechleiter@uthscsa.edu

-----  
Name                      Badge#                      Contact#                      Email

## Authorized Users

While an important part of your SOP, this list of all Authorized Users, including their badge number, should appear as Appendix C of this document and users should sign the document to verify that they have read and understand the Standard Operating Procedure for the LCA. Appendix C can be updated as frequently as needed.

## Incidental Personnel

In addition to the authorized primary users, incidental personnel may be in the room at the time of the experiment. These personnel are not trained on the system, nor have they gone through the laser safety program. They will be observing only and not involved in the operation of the system. They will not be in the room during any alignment or maintenance of the lasers.

## Normal Laser Operation

Nikon Sweptfield confocal contains the following lasers:

Class	Type	Make	S/N	Wavelength (nm)	Power output (mW)
3B	Solid state	Agilent	N1245AL14	405	50
3B	Solid state	Agilent	N1245AL24	488	50
3B	Solid state	Agilent	N1245AL34	561	50
3B	Solid state	Agilent	N1245AL44	640	50

All lasers and their power supplies are integrated into a wheeled unit that is placed on the floor under the anti-vibration table. Laser output is delivered to the Nikon Ti microscope via a fiber optic cable. The S/N for the integrated laser launch unit is N1245B ATO-3598 SER US51360103.

All lasers are used as the excitation sources for confocal microscopy.

## Eyewear section

Wavelength specific eyewear will be used by field service engineers during alignment.

## Alignment Hazard Control

All lasers are aligned by service engineers and are not adjusted by users or facility staff.

## Laser Hazard Control

1. Access to laser control room is restricted to trained personnel. The room is secured by a card-access door reader. All access is approved by the campus police.

2. The entire system is enclosed by laser proof curtain.
3. Operation of all Class 3B lasers is only accessible through the system software. Computer access is restricted to trained users with unique login names and passwords.
4. A "Laser in Use" warning sign is mounted above the entrance of the enclosed system booth.
5. All lasers are enclosed with no open beam throughout the path.
6. All users are trained in the operation of the NSTORM and the proper use and care for the lasers.

Primary users are trained in the operation of the Sweptfield confocal and the proper use care for the lasers.

### **Control of Additional LCA hazards**

Indicate other possible hazards associated with the lasers in your LCA.

### **Associated Chemical Hazard Control**

List chemicals used in this LCA include a list MSDS numbers, or attach MSDS's to the end of the document (not required to attach MSDS sheets). If you prefer, provide the chemical list section from your Project Review Document as an Appendix. Indicate in this section if there is a registered Satellite Waste Accumulation Area in the LCA and where it is located. Discuss any site-specific chemical hazards for this LCA in this section.


### **Emergency Procedures**

1. Shut down the laser system. Use "emergency stop" button if equipped.
2. Provide for the safety of the personnel, i.e. first aid, CPR, etc.
3. If a fire has been created as a result of the laser, follow appropriate procedures to put out the fire:
  - Within the surgical field, douse with sterile water
  - For fires in other areas, utilize an appropriate class fire extinguisher to extinguish the fire
4. Obtain medical assistance. In the event of a suspected eye injury, have the injured person keep their head upright and still to restrict any bleeding in the eye.
  - For life-threatening injuries (major burns, cardiac arrest following electrocution), dial 911 for immediate medical assistance.
  - For non-life threatening injuries (laser eye injuries, minor skin burns), employees should be evaluated by a physician as soon as possible. **Do not allow anyone with a potential laser eye injury to drive themselves.**

Emergency Dispatch

911 (from campus phone)





This SOP has been adapted with permission from Argonne National Laboratories.

**User Name**

Aruna Jaiswal  
Benjamin Jordan  
Hema Manogna Gudlavalleti  
Jessica Johnson  
Juli Bai  
Lauren Schneider  
Linda Phillips  
Manuel Riquelme  
Nozomi Tomimatsu  
Sadiya Ahmad  
Sri Balasubas Muralimanoharan  
Suman Kanji  
Victoria Alers  
Yan Zeng  
Yilun Deng

**User Login Email**

jaiswala@uthscsa.edu  
jordanb1@uthscsa.edu  
gudlavalleti@uthscsa.edu  
johnsonj20@livemail.uthscsa.edu  
baij@uthscsa.edu  
schneiderl1@uthscsa.edu  
phillipsl@uthscsa.edu  
riquelme@uthscsa.edu  
tomimatsu@uthscsa.edu  
ahmadsn@livemail.uthscsa.edu  
muralimanoha@uthscsa.edu  
kanji@uthscsa.edu  
alers@uthscsa.edu  
zengy1@uthscsa.edu  
dengy@livemail.uthscsa.edu